

# Technical Specifications

MolecuLight *i:X*<sup>™</sup>

---



---

## Capacity

32GB

---

## Display Screen

4-inch (diagonal) widescreen display with Multi-Touch IPS technology  
1136-by-640-pixel resolution at 326 pixels per inch  
800:1 contrast ratio (typical)  
500 cd/m<sup>2</sup> max brightness (typical)  
Fingerprint-resistant oleophobic coating

---

## Size & Weight

**Height:** 21 cm  
**Width:** 9.5 cm  
**Depth:** 4 cm  
**Weight:** 398 g

---

## Camera & Video

5-megapixel camera  
Autofocus  
Tap to focus video or still images  
*f*/2.4 aperture  
Video stabilization  
Image format captured: JPG  
Video format captured: MOV (1080p)

---

## Violet Wavelength LEDs

Emits safe violet wavelength light (405 nm)

# Technical Specifications continued

---

<b>Fluorescence Emission Filters</b>	500 - 545 nm 600 - 665 nm
<b>Sensors &amp; Indicators</b>	Camera Sensor Ambient Light Status LED Range Finder LED System Status LED Display Screen Battery Indicator
<b>External Buttons &amp; Connectors</b>	Power Button Display Screen Power Button Home Button Display Screen Charging Port MolecuLight <i>i:X</i> Charging Port Rocker Switch (Fluorescence Imaging Mode™ and Standard Imaging Mode™)
<b>Power &amp; Battery</b>	Built-in rechargeable lithium-ion battery MolecuLight Power Cable used to charge MolecuLight <i>i:X</i> Imaging Device and Display Screen
<b>Operating System</b>	iOS 9.3.5
<b>Built-in Apps</b>	iX Imaging App (Version 1.1)
<b>Operating Conditions</b>	<b>Temperature:</b> 10 °C - 30 °C <b>Relative Humidity:</b> 10 - 70 % <b>Atmospheric Pressure:</b> 89 - 102 kPa
<b>Cleaning</b>	Compatible with disinfecting wipes containing 70% isopropyl alcohol

---

For more information about the MolecuLight *i:X*, please refer to the MolecuLight *i:X* User Manual.



©2021 MolecuLight® Inc. All Rights Reserved. PN 1395 Rev 1.2

The MolecuLight *i:X*™ Imaging Device is approved by Health Canada and has CE marking for sale in Canada and the European Union. The MolecuLight *i:X* Imaging Device has received FDA clearance.

MolecuLight® is a Registered Trademark in Canada, the US, and the EU.